

For life science research only
Not for use in diagnostic procedures



iMatrix-511

Product No. 892 011 350 µg
Product No. 892 012 1,050 µg

Version 005
Store at 2-15 °C
Protect from light

Background Information

Laminin-511 is well known to bind to the integrin $\alpha 6 \beta 1$ which is located on the cell surface. iMatrix-511 is recombinant Laminin511-E8 fragments.

Content

Recombinant Human Laminin511-E8 Fragments

Amount

175 µg / tube (892 011: 2 tubes, 892 012: 6 tubes)

Concentration

0.5 mg / mL

Form

Liquid solution (solvent: PBS(-))

Product Information

iMatrix-511 is recombinant human Laminin511-E8 fragments expressed by CHO-S cell (Life Technologies).

Storage and Stability

The liquid solution is stable at +2 to +15 °C until the expiration date printed on the label.

Protect from light.

iMatrix-511 is stable at 4 °C

for 2 years from the manufacturing date.

Activity

The dissociation constant of the binding activity with integrin $\alpha 6 \beta 1$ is under 10 nM.

Application

iMatrix-511 is able to use as cell culture substrate for various cell types including ES/iPS cells.

Procedure

1) Dilute the solution with sterile PBS(-). Coat dishes with 0.5 µg/cm².

* For example, for one well of a 6-well plate (9.6 cm² /well), add 9.6 µL of iMatrix-511 (4.8 µg) in 1.99 mL of PBS(-).

Add 2 mL of diluted iMatrix-511 solution to the well.

2) Incubate for 1 h at 37 °C, 3 h at room temperature, or over night at 4 °C.

3) Remove remaining fluid from the coated surface. No rinse is needed.

4) Immediately plate the cells at desired density.

* Don't allow the plate to dry.

* Briefly spin down all liquid in the tube before use.

* Avoid repeated freeze-thaw cycles.

References

Ido H et al. The requirement of the glutamic acid residue at the third position from the carboxyl termini of the laminin gamma chains in integrin binding by laminins. *J. Biol. Chem.* **282** (15): 11144-54, 2007

Taniguchi Y et al. The C-terminal region of laminin beta chains modulates the integrin binding affinities of laminins. *J. Biol. Chem.* **284** (12): 7820-31, 2009

Miyazaki T et al. Laminin E8 fragments support efficient adhesion and expansion of dissociated human pluripotent stem cells. *Nat. Commun.* **3**: 1236, 2012

Nakagawa M et al. A novel efficient feeder-free culture system for the derivation of human induced pluripotent stem cells. *Sci Rep.* **4**: 3594, 2014

Doi D et al. Isolation of Human Induced Pluripotent Stem Cell-Derived Dopaminergic Progenitors by Cell Sorting for Successful Transplantation. *Stem Cell Reports.* **2** (3): 337-50, 2014

Takashima Y et al. Resetting Transcription Factor Control Circuitry toward Ground-State Pluripotency in Human. *Cell.* **158** (6): 1254-69, 2014

Fukuta M et al. Derivation of mesenchymal stromal cells from pluripotent stem cells through a neural crest lineage using small molecule compounds with defined media. *PLoS One.* **9** (12): e112291, 2014

Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

Contact and Support

E-mail: protein-info@nippi-inc.co.jp

Research & Development Project Manager, NippiResearch
Institute of Biomatrix, Protein Engineering Office, Nippi,
incorporated
1-1-1 Senjumidori-cho, Adachi-ku, Tokyo 120-8601, JAPAN