

## Cell lysis for protein extraction\*

- 1. Carefully remove decant culture medium from adherent wells.
- 2. Wash all cells twice with cold PBS.
- 3. Add ice-cold RIPA buffer (with inhibitors) to the cell wells. Use  $\sim$ 600  $\mu$ L of buffer per well (12 well plates), it is important that the fibers are covered in buffer.
- 4. Keep on ice for 20 minutes, swirling the plate occasionally for uniform spreading.
- 5. Collect the lysate and transfer to a cold microcentrifuge tube. Centrifuge samples at 4°C at 14 000 x g for 15 minutes to pellet the cell debris.

  Note: To increase yield, sonicate the pellet for 30 seconds with 50% pulse.
- 6. Gently remove the tubes from the centrifuge and place on ice, aspirate the supernatant and place in a fresh tube kept on ice.

<sup>\*</sup>Suggested procedure, please adjust according to your experimental needs.