

Detachment Method using EDTA for human iPS cells cultured on iMatrix[™]-511 (laminin-511E8)

(In the case of 6 well plate)

- 0) Culture human iPS cells for 7 days on the plate pre-coated with 0.5 μg/cm² iMatrixTM-511
- 1) Aspirate the used medium.
- 2) Wash the cells with 2 ml/well of 5 mM EDTA/PBS(-) twice.
- 3) Incubate the cells with 1 ml/well of 5 mM EDTA/PBS(-) at 37°C for 15 min
- 4) Aspirate EDTA/PBS(-).
- 5) Add 1 ml/well of StemFit medium containing Y27632. Pipet repeatedly to detach and dissociate the cells into single cells.
- 6) Count the number of the cells and seed the cells accordingly.

^{*} This protocol is NOT the final version. We may possibly modify the protocol later.

^{**} The cells maintained high viability, high proliferation rate, undifferentiated state and normal karyotype after 10 passages with this protocol.

^{***} Incubation time depends on the cell line. Incubation time may be shortened to 10 minutes.

^{****} Protocol should work with ES Cells but currently untested.