

Artificial serum

Serum substitute

(Supplement for serum-free culture media)



Culture media
supplement

Artificial serum
Animal-free




- For human Mesenchymal Stem Cells, human Fibroblasts, human Lymphocytes and other cell lines.
- A small amount of Artificial serum, 5% or less, supports cell growth just as well as FBS.
- No lot-to-lot variation

Artificial serum (Animal-free)



- Artificial serum is developed as a substitute of serum (e.g. FBS).
- Additional "Artificial serum" to serum-free media enhances cell growth.
- Recommended concentration is 5% or less. Optimal concentration should be determined for specific use and cell type.
- The components of Artificial serum are animal-free.



1. Adhesion cells (Fibroblast)



2. Suspension cells (DAUDI)



Code (CSTI / NIPRO)	Product name	Description	Storage	Volume
A2G20P2CC / 87-082	Artificial serum (Animal-free)	Animal-free Supplement	below -20°C	20mL (P)

(P) PET bottle

※ For Research use only.



1-16-16 Seikaen, Aoba-ku, Sendai,
989-3162, Japan
URL: <http://www.cstimedia.com>
E-mail: csti-info@ml.nipro.co.jp



NIPRO CORPORATION
3-26, Senriokashinmachi, Settsu,
Osaka, 566-8510, Japan

Official Distributor for North America

Iwai North America Inc.



Tel: 650-486-1541
orders@iwai-chem.com
www.iwaichem.com