


BMPro[®]-W

Serum free medium for Antibody/Recombinant protein production



Suspended CHO Cells

- Hydrolysates free
- Serum component free
- Low protein contents


BMPro[®]-W

BMPro-W is serum-free complete medium formulated for CHO cells thriving in suspension culture systems. Tailored for optimal antibody and recombinant protein production, this medium is devoid of serum, hydrolysates, and maintains a low protein content. To enhance performance, simply supplement with 8mM L-glutamine.




Proliferation of CHO cell in suspension culture

Fig.1 Cell proliferation curve



Seeding density : 2x10⁵ cells/mL

Fig.2 Viability



Catalog# (CSTI / NIPRO)

Product Name

Description

Size (Packaging)

2502P10 / 87-015

BMPro-W

Serum free media for antibody and recombinant protein production (Optimized for CHO cell in suspension culture system) L-glutamine free

1000mL (P)

Supplement

Catalog# (CSTI / NIPRO)

Product Name

Description

Size (Packaging)

1104 / 87-885

200mM
L-Glutamine solution

200 mM L-glutamine Solution

100mL (P)

(P) PET bottle

※ This product is for research use only. It cannot be used for therapeutic or diagnostic purposes.



株式会社細胞科学研究所
仙台市青葉区西花苑1丁目16番16号
URL: <http://www.cstimedia.com>
E-mail: csti-info@ml.nipro.co.jp



ニプロ株式会社
大阪市北区本庄西3丁目9番3号

Official Distributor for North America

Iwai North America Inc.



Tel: 650-486-1541
orders@iwai-chem.com
www.iwaichem.com