



Live transport system

Dish Cover Compatibility

Updated: Jan. 30, 2024

	φ35 dish		φ60 dish		φ90 ~100 dish	
Starting volume (mL)	4~6		14~17		44~55	
Make	Fit	Model	Fit	Model	Fit	Model
SUMILON	○	MS-11350			○	MS-13900
IWAKI	○	1000-035 3000-035 3971-035 4971-011	○	3010-030 1010-060	○	1020-100 3020-100
FALCON	×	353001	○	353002	○	430293
NUNC	○	174913				
CORNING	○	430165	○	430166	○	430591 430293
BD Falcon			○	353002	○	
greiner	○	627160	○	628160	○	664160-013
NEST	○		○		○	
STEM						
SARSTEDT	○	831 800 002	○	83.3901.002	○	83.3902.300
	×	83.3900				
Matsunami Glass Ind., Ltd.	○					
VIOLAMO	○	2-8590-01	○	2-8590-02	○	2-8590-03
CellSeed	○	UpCell	○	UpCell HydroCell	○	UpCell HydroCell
Legend: ○= tested and compatible ×= tested and not compatible Blank cell: not tested						