



Live transport system

Microplate Cover Compatibility

Updated: Jan. 30, 2024

	6 Well Plate		12 Well Plate		24 Well Plate		96 Well Plate	
Starting volume per well (mL)	9		4		2		0.2	
Make	Fit	Model	Fit	Model	Fit	Model	Fit	Model
CellSeed	○	UpCell®	×	UpCell®	○	UpCell HydroCell	○	UpCell®
	×	HydroCell						
Corning®/Elplasia®					×	4441		
Corning®/Costar®	×	3516	×	3513	×	3526	×	3596
							○	3474
eppendorf	○	0030 720.016	○	0030 721.012	○	0030 722.019	○	0030 730.011
FALCON®	○	353046	○	353043	○	353047	○	351172 353072
Greiner	○	657160	○	665180	○	662160		
IWAKI	○	3810-006	○	3815-012	○	3820-024	○	3860-096
					○	3820-024 (Glass bottom)		
MatTek					○	P24G-1.5-13-F		
NEST	×		×		×		○	
NUNC™	○	140675	×	150628 174910	○	142475	○	167008 174907 174908
Sarstedt	×	83.3920	×	83.3921	×	83.3922	○	83.3924
STEM	×	Flat Bottom			×	Flat Bottom	○	Flat, Round & V bottoms
SUMILON®	○	MS-80060	○	MS-80120	×	MS-80240	○	MS-8096F
VIOLAMO	×	2-8588-01	×	2-8588-02	×	2-8588-03	○	2-8588-05

Legend: ○ = tested and compatible × = tested and not compatible Blank cell: not tested